

MAGNETIC SORTING OF TRANSFECTED CELLS USING THE MACSELECT  
SYSTEM

**Materials:**

MACSelect K<sup>k</sup>.II plasmid

OctoMACS magnet

MS separation columns

PBE:

PBS supplemented with 0.5% BSA and 5 mM EDTA

PBE must be degassed before use

**Procedure:**

1. Transfect cells by method of choice on 100 mm plate.

Cotransfect your plasmid of interest along with the MACSelect K<sup>k</sup>.II plasmid (encodes mouse MHC class I H-2K<sup>k</sup> protein with truncated cytoplasmic domain to avoid signal transduction). Use a 1:1 ratio.

OR

Clone your gene of interest into the MACSelect K<sup>k</sup>.II plasmid and transfect into your cell line.

2. Incubate 24-48 hours (determine optimal time for your cell type- best expression occurs after one cell doubling).
3. Harvest less than or equal to  $1 \times 10^7$  transfected cells (out of a maximum number of  $2 \times 10^8$  total cells):
  - a. Wash cells with PBS, removing all dead cells.
  - b. Add 500  $\mu$ l trypsin and allow to trypsinize in incubator.
  - c. Stop trypsinization by adding 100  $\mu$ l of 100% serum.  
To avoid clumping of HeLa cells: Add 5 ml of PBE to trypsinized cells and resuspend. Check under the microscope that there are no clumps. Spin down the cells and remove the supernatant. Add 600  $\mu$ l PBE and resuspend.
  - d. Add 80  $\mu$ l of MACSelect K<sup>k</sup> Microbeads (binds to cells expressing the MACSelect K<sup>k</sup>.II- encoded protein).
  - e. Rock plate gently by hand to disperse. Incubate 15 min at RT, rocking twice more during the incubation.
  - f. Add PBE to adjust the final volume to 2 mls (add 1320  $\mu$ l) and resuspend cells completely. Verify cells are not clumped. If they are, separate by passing several times through an 18 gauge needle or a 1 mL pipette tip.

4. Magnetic separation:

- a. Attach OctoMACS Separation Unit to the MACS MultiStand.
- b. Put MS column on the OctoMACS Separation Unit. Place a waste tube under the column.
- c. Apply 500  $\mu$ l PBE to column and let run through into collection tube.
- d. Save 100  $\mu$ l aliquot of cells (for “pre-column” control), and apply the rest to the column 500  $\mu$ l at a time. Resuspend each aliquot well before applying to the column.
- e. Wash the column 4 times with 500  $\mu$ l PBE.
- f. Place column in collection tube.
- g. Add 1 ml PBE, and flush out cells with plunger into collection tube.

5. Fluorescence detection:

- a. Take 100  $\mu$ l of cells in PBE (from the “pre-column” and the “post-column” samples) and add 10  $\mu$ l of H-2K<sup>k</sup> FITC.
- b. Incubate in the dark at 4C.
- c. Wash with 1 ml PBE, spin, and remove supernatant.
- d. Resuspend pellet without additional buffer.
- e. Transfer to glass slide and visualize by microscopy.
- f. Determine percent enrichment.
- g. Cells can be enriched again on a fresh column to give higher percent enrichment if necessary.

**Reference:**

Miltenyi Biotech, [www.miltenyibiotec.com](http://www.miltenyibiotec.com)

**Submitted by:** Carina Holmberg and Sandy Westerhide