

## CAT ASSAY PROTOCOL

### PROTOCOL #1

1. Wash plates 2x in PBS.
2. Scrape into microfuge tube in 1 ml PBS.
3. Spin out cells and aspirate off PBS.
4. Resuspend in 100  $\mu$ l 0.25 M Tris pH 8. Vortex.
5. Freeze in dry ice/EtOH. *Can be stored at this point.*
6. When ready, freeze/thaw in 37°C (until thawed) to -70°C ice-bath, back and forth 3 times. On third time, freeze again at -70°C.
7. Let thaw at room temp.
8. Spin in cold microfuge 5 min. Leave pellet and extract in tube.
9. For protein determination:  
For each tube mix:  
798  $\mu$ l 50 mM tris pH 8  
2  $\mu$ l extract  
200  $\mu$ l Bradford dye  
  
Vortex.  
Read O.D. at 595 nm vs. blank (as above with no extract).
10. Read concentrations off standard curve.  
On Gregg Williams' standard curve, values are per  $\mu$ l; since we use 2  $\mu$ l extract, divide number in half to get  $\mu$ g/ml.
11. Calculate amount needed for 100  $\mu$ g and amount of 0.25 M tris pH 8 to bring total volume to 149  $\mu$ l.
12. For each reaction:  
149  $\mu$ l extract in 0.25 M Tris 8.0  
1  $\mu$ l <sup>14</sup>C chloramphenicol  
20  $\mu$ l 4 mM Acetyl CoA T--6.6436 mg in 2 mls  
50 mM Tris pH 8 (or 3.3218 mg in 1ml)  
  
(MW = 809.6 + 20.85 (Li) = 830.49)
- Make rxn mix of labeled chloramphenicol & acetyl CoA for n + 1; add 21  $\mu$ l to each tube containing the Tris, then add extract last to start the rxn (on ice until start). Incubate at 37°C for 1 hour.
13. Add 1 ml cold ethyl acetate to each; vortex.
14. Spin 5 seconds in cold.
15. Take off top layer carefully (completely avoid interphase) into new tube.
16. Spin vac to dryness, ~1/2 hour.
17. For chromatography tank, mix in a beaker:  
95 ml chloroform  
5 ml methanol  
  
Carefully place in the bottom of the tank with a 25 ml pipette - do not splash onto sides of tank. Close tightly and let equilibrate while samples are drying down.
18. Resuspend pellets in 20  $\mu$ l ethyl acetate and spot onto TLC sheet one inch from bottom. Place in tank and allow to chromatograph until liquid is ~2" from top.
19. Remove carefully, let dry to completeness in hood.
20. Place on film at room temperature, no screen.

## PROTOCOL #2

1. Transfect cells with 5 µg plasmid DNA encoding activator, 5 µg reporter plasmid containing CAT gene, 2 µg internal control plasmid for transfection efficiency (e.g. plasmid encoding luciferase gene or B-galactosidase gene) and 8 µg of carrier DNA.
2. Harvest cells within 48 hrs after transfection and freeze cell pellets quickly in dry ice. Then store the cell pellets at -80°C overnight or lyse the cells by adding 4 volumes of buffer C to cell pellets, put on ice for 10 min, spin 13,000 rpm 15 min, move the supernatant to a new tube and freeze rapidly in dry ice/methanol. Store at -80°C or do a BioRad protein assay directly.
3. Combine 1 µl 25 µCi/ml [<sup>14</sup>C] chloramphenicol, 20 µl 4 mM acetyl CoA and 32.5 µl 1M Tris.Cl, pH 7.4. Add 12.5 µg cell extract and add ddH<sub>2</sub>O to 150 µl, mix and incubate at 37°C for 1 hr.
4. After incubation add 1 ml ethyl acetate to the reaction and vortex. Spin 1 min at 4°C, transfer top layer to a new tube, dry in Speed-vac for 30 min, and resuspend in 15 ml ethyl acetate.
5. Spot samples evenly on a TLC plate and develop in a chromatography tank containing 19:1 of chloroform : methanol. The tank should be equilibrated 2 hrs prior to developing. Run the sample until the solvent is close to the top of the plate.
6. Remove TLC plate from the tank, air dry and expose to X-ray film at room temperature. No intensifying screen is needed.

**Reference:** Gorman, C.M., Moffat, L.F., and Howard, B.H., 1982, MCB 2: 1044-1051