

LUCIFERASE ASSAY: (PROMEGA SYSTEM)
(24 Well Plate)

Cell Lysis

1. Thaw 5X PLB and prepare appropriate amount of 1X PLB (100 μ l per well) by diluting in dH₂O
2. Remove media from cells by aspiration
3. GENTLY wash with 1X PBS (perhaps should be skipped for 293T cells)
4. Add 100 μ l 1X PLB to each well. Place on shaker at RT for 30 min.

Reagent Preparation

5. Prepare Luciferase assay reagent II
Resuspend lyophilized Luciferase Assay Substrate in 10 mL of Luciferase Assay Buffer II and label "LARII"
(stable for one year at -80)
* need 100 μ l per assay, plus about 2 mLs to tube volume for luminometer

Luminescence measurement

6. Read samples in Molecular Devices "LMax" luminometer
Dispense 20 μ L cell lysate into white opaque luminometer microtiter plate
Prime P injector (7 injections) with LARII
Inject 100 μ l LARII per well with P-injector.
Read with a 2 sec premeasurement delay, 10 sec integration to detect firefly luciferase.

Wash injector

Place P injector tube into a 50ml Falcon tube containing ddH₂O. Select **Wash Injectors...** from the **Control** menu. Accept the default settings (30 injections) and click **OK**.

Repeat this wash step with 70%EtOH (30 injections). Perform a final "Wash" that will just dry the injector tube.

Protein quantitation

Quantitate protein concentrations for each well by Bio-Rad assay.

Data analysis

Normalize luciferase numbers to protein concentration numbers.

Submitted by: Sandy Westerheide