

FREEZING CELLS

Materials:

Fetal Bovine

Sera

DMSO (Sigma D-5879)

1.8ml cryovials (Nunc/Nalge 368 632)

1xPBS

media/10%FBS

1xtrypsin

Freezing Solution

90% FBS

10% DMSO

Procedure:

1) Grow cells until ~80% confluency.

80% confluency	amount of freezing solution	number of cryovials
25 cc flask	1 ml	2
100mm plate	1 ml	2
150mm plate	2ml	4

- 2) Wash cells with 1xPBS, trypsinize (if necessary), then add 10ml fresh media/10%FBS and transfer to a 15ml conical tube.
- 3) Spin cells 2000rpm 2 min.
- 4) Label cryovial with:
 - cell line name
 - date
 - your initials
- 5) Discard media and resuspend cell pellet in freezing solution (mix well).
- 6) Aliquot 0.5ml into the labeled cryovials.
- 7) Cells need to freeze slowly. Place vials in the top chamber rack of the blue liquid nitrogen tank **overnight** to freeze. Vials will be above the level of the liquid nitrogen and will freeze gradually.
- 8) Put the frozen vials into a storage box in one of the tanks.
- 9) Fill in the information in the logbook and database record sheet.