

### MULTIPRIME DNA LABELLING

1. Follow the protocol specified by Amersham, Inc.
2. Circular DNA is linearized by restriction digestion then heated to inactivate restriction enzyme; DNA is diluted to 10 ng/ $\mu$ l in TE.
3. Linearized or gel-isolated fragment DNA is denatured by heating to 95-100°C for 2 minutes then quickly chill on ice.
4. Use 2.5  $\mu$ l of DNA (25 ng) per reaction.
5. Add 15  $\mu$ l multiprime buffer solution 1 (supplied by Amersham and containing random primers, nucleotides and buffer).
6. Add H<sub>2</sub>O to 43  $\mu$ l, then 5  $\mu$ l of [<sup>32</sup>P]dCTP (>2000 Ci/mmol).
7. Add 2  $\mu$ l of multiprime solution 2 (containing Klenow fragment of *E. coli* DNA polymerase I). Mix by gently tapping tube.
8. Place at room temperature for 3-16 hours.
9. Check incorporation by TCA precipitation of 1  $\mu$ l and store at -20°C.
10. Because incorporation is 90%, no separation of unincorporated dNTPs is necessary.