

ELECTROPHORESIS OF RNA THROUGH GELS CONTAINING FORMALDEHYDE

1. Prepare gel-running buffer and formaldehyde.

a. 5 X Gel-running buffer:

0.2 M morpholinopropanesulfonic acid (MOPS) (pH 7.0)

50 mM sodium acetate

5 mM EDTA (pH 8.0)

This buffer yellows with age if exposed to light or autoclaved. Discoloration does not affect its performance appreciably.

b. Formaldehyde (F.W.=30.03) is usually obtained as a 37% solution in water (12.3 M). Check that the pH of the concentrated solution is greater than 4.0. The concentrated solution should be handled and stored in a chemical hood.

2. Prepare the gel by melting agarose in water, cooling to 60°C, and adding 5X gel buffer and formaldehyde to give 1X and 2.2 M final concentrations, respectively. (One part of stock formaldehyde solution should be diluted with 4.6 parts of agarose solution).

*Note: The fractionation properties of formaldehyde gels with different agarose concentrations have been determined by Lehrach et al. (1977).*

3. Prepare the sample by mixing the following in a sterile Eppendorf tube:

RNA (up to 20 µg)	4.5 µl
5 X gel-running buffer	2.0 µl
formaldehyde	3.5 µl
formamide	10.0 µl

Incubate at 55°C for 15 minutes.

*Note: Formamide oxidizes readily in air and should be deionized by passage through a mixed-bed resin (Bio-Rad AG 501-X8) until its pH is neutral. It is then recrystallized at 0°C and stored at -20°C in small aliquots in tightly capped tubes.*

4. Add 2 µl of sterile loading buffer.

Loading buffer:

50% glycerol

1 mM EDTA

0.4% bromophenol blue

0.4% xylene cyanol

5. Load the RNA samples onto the gel. Restriction fragments of DNA are

convenient molecular-weight markers. They should be treated and run exactly as the RNA samples. Labeled DNA markers or markers that will be detected by a labeled hybridization probe are preferred because the ethidium fluorescence of formaldehyde-denatured nucleic acids is weak. If such markers are impossible to obtain, the following protocol may be used:

- a. Apply to the gel sufficient DNA to give at least 50-100 ng per band.
- b. Cut the lanes containing the markers from the gel before the alkaline hydrolysis step (see below), and wash with four or five changes of water for 2 hours.
- c. Wash with two changes of 0.1 M ammonium acetate for 1 hour.
- d. Stain for 1 hour with 0.5  $\mu\text{g/ml}$  of ethidium bromide in 0.1 M ammonium acetate, and 0.1 M  $\beta$ -mercaptoethanol.
- e. Destain for 45 minutes with a solution of 0.1 M ammonium acetate and 0.01 M  $\beta$ -mercaptoethanol.

#### Transfer of Formaldehyde-denatured RNA to Nitrocellulose

1. After electrophoresis is complete, soak the gel for 5 minutes in several changes of water.

*Note: Gels containing formaldehyde are less rigid than nondenaturing agarose gels. Care must be exercised in handling them.*

2. Soak the gel in an excess of 50 mM NaOH and 10 mM NaCl for 45 minutes at room temperature.

*Note: The partial alkaline hydrolysis improves the transfer of high-molecular-weight RNA.*

3. Neutralize the gel by soaking for 45 minutes at room temperature in 0.1 M Tris Cl (pH 7.5).
4. Soak the gel for 1 hour in 20X SSC.
5. Transfer the RNA to nitrocellulose by the method described for Southern blots. The transfer is complete in 3-4 hours.
6. After transfer is complete, wash the filter in 3X SSC, dry in air for 1-2 hours, and bake for 3-4 hours at 80°C under vacuum.

#### Reference:

Procedure from Maniatis cloning manual.

*Lehrach, H., D. Diamond, J.M. Wozney and H. Boedtker. 1977. Biochemistry 16, 4743.*