

CALCIUM PHOSPHATE MEDIATED TRANSFECTION

Materials:

- | | | |
|----|--|----------------|
| 1) | 2 M CaCl ₂ .2H ₂ O | 29.4 g/100mls |
| | | Filter |
| 2) | 100X PO ₄ | |
| | 70 mM Na ₂ HPO ₄ | 0.994 g/100mls |
| | 70 mM NaH ₂ PO ₄ | 0.966 g/100mls |
| | | Filter |
| 3) | 2X HBS | |
| | 10 g/liter Hepes | |
| | 16 g/liter NaCl | |
| | make 100 mls, pH to 7.15, Filter. | |

Procedure: (Amounts refer to a transfection of one p100 of cells)

- 1) Precipitate 10-20 µg of DNA. DNA must be extremely pure. 2 CsCl gradients are highly recommended.
- 2) In hood, remove ethanol.
- 3) Without disturbing pellet, fill tube with cold 100% ethanol to sterilize. Remove 100% ethanol carefully. Rinse cap.
- 4) Invert to dry. Drying occurs rapidly in hood. Do not overdry.
- 5) Resuspend pellet in 440 µl H₂O.
- 6) Add 60 µl 2 M CaCl₂.
- 7) Incubate at 37°C >4 hrs to ensure that the DNA is fully dissolved. This DNA/calcium solution is called solution B.
- 8) Make solution A (in hood).
500 µl 2X HBS
10 µl 100X PO₄.
- 9) Very slowly, dropwise, add solution B to solution A.
- 10) A very fine, cloudy precipitate should form. Ideally no particles should be visible, only a grayish haze. If extensive clumping occurs, start over. Clumping can be due to doing step 9 too fast or by failing to completely resuspend the DNA in step 7. Let precipitate sit in hood for 15-30 minutes to allow precipitate formation to go to completion.
- 11) Pipette precipitate vigorously to break up any minor clumping that might have occurred.

Steps 12 -- differ for 293 and Hela Cells.

293 cells

- 12) Wash cells (~85% confluent) 2X with warm DME. (Do this while ppt is sitting)
- 13) Add 7mls DME to plate.
- 14) Dropwise, add ppt to plate. Tilt plate to mix.

- 15) Let incubate O/N, 37°C, 5% CO₂.
- 16) Remove media.
Replace with 8mls normal growth media (containing serum).
- 17) For transients harvest 24-48 hours later.
For stable lines, split 1:5 to 1:10 24 hours later. Incubate 24 hours and then select. (400 µg/ml G418).

Hela cells

- 12) Wash cells (~85% confluent) 2X with warm DME. (Do this while ppt is sitting)
- 13) Add 1 ml ppt directly to cells. (No DME)
- 14) Incubate 30', 37°C, 5% CO₂.
- 15) Add 7 ml DME + 2% Fetal Calf Serum.
- 16) Incubate 6-8 hours, 37°C, 5% CO₂.
- 17) Remove media.
Shock with 2ml DME + 25% glycerol.
Shock for 2-3 minutes.
- 18) Wash with warm DME 2X.
- 19) Add 8ml normal growth media (containing serum).
- 20) For transients, harvest 24-48 hours later. For stable lines, split 1:2 to 1:5 24 hours later.
Incubate 24 hours, then add selective media (400 µg/ml G418).

Comments:

Pooling of 293 cells can be done after 2-3 weeks. Helas, >1 month. While selecting change media when dead cells become the majority or if the media becomes orange.