

## DRUG SCREEN USING MAMMALIAN CELLS

This procedure has been developed for large-scale drug screens (up to 6000 compounds) using mammalian cells. We used HeLa HSE-Luc cells to detect for the induction of the heat shock response. We pool the drug compounds into a “master mix” before adding to the cells (6 drugs combined at 5 $\mu$ M each). Thus, 6000 drugs can be initially screened using eleven 96-well plates of cells. We then analyze for induction of luciferase expression using Promega’s Bright-Glo Luciferase Assay detection system.

### Materials:

expression cells (containing reporter construct)  
hemacytometer  
8-well multi-channel pipettor  
sterile multi-channel reservoir trays  
96-well tissue culture plate (Corning #3596 flat-bottom cell culture)  
Round bottom 96-well plate for drug “master plate” (Falcon #3076 U-bottom)  
96-well plate for luminometer (Corning #3912 white polystyrene flat bottom)  
Luminometer (we have a Molecular Devices *Lmax* luminometer)

### Reagents:

cell culture media, 1xPBS, 1xTrypsin, G418 (200 $\mu$ g/ $\mu$ l)  
DMSO (Sigma D-5879)  
Bright-Glo Luciferase Assay reagent (Promega E2620)  
Drugs (we obtained the drugs as 100 $\mu$ g lyophilized on 96-well plates)  
Celesterol (100 $\mu$ M) for positive control  
CdCl<sub>2</sub> (1mM) for positive control

To split cells to 96-well plate:


1. Grow cells in tissue culture flask or plate (HeLa HSE-Luc cells grow in DMEM/10%FBS with G418 geneticin at 200 $\mu$ g/ml). Allow to grow until ~80% confluent.
2. Split cells to tissue culture 96-well plate as indicated below when ~80% confluent (and ~24 hours before drug treatment):
  - 1) remove media
  - 2) rinse plate with 1xPBS
  - 3) add 1xTrypsin to cover the plate
  - 4) incubate at 37°C ~5min
  - 5) add 8ml media to plate; collect cells and transfer to a 15ml conical tube
  - 6) centrifuge to pellet cells 2000rpm for 2min
  - 7) resuspend pellet in volumes suggested below and mix
  - 8) remove 10 $\mu$ l to count cells using the hemacytometer.

Hela Hsp70-Luc cells if ~80% confluent		
plate size	resuspension volume	approx. number of cells
100mm		
150mm	1.5-2.0 ml	4,500,000 cells

- 9) For one 96-well plate you will need 700,000 cells in 10ml (to have a little extra).  
for example, if 200,000 cells/ml:  
 $200,000 \text{ cells/ml} / 700,000 \text{ cells} = 0.286 \text{ ml cell suspension} + 9.71 \text{ ml media} = 700,000 \text{ cells/10ml}$
- 10) In a tissue culture hood, combine cells and media in a sterile multi-channel pipet reservoir tray and mix well.
- 11) plate 100 $\mu$ l per well to a 96-well plate (7000cells). *Note: you will need extra wells of cells for controls. We did controls in duplicate on the same plate.*
- 12) cells should be ~80% confluent after 18-24 hours.

### Drug Treatment

1. Resuspend lyophilized (100 $\mu$ g) drugs in 100 $\mu$ l DMSO (to bring the concentration to 1 $\mu$ M/ $\mu$ l). Dissolve the drugs by pipetting up-and-down several times. Once resuspended the plates may be stored at -20 $^{\circ}$ C (the DMSO will take ~1 hour to thaw at room temp once frozen).
2. Prepare a master mix plate using a sterile round-bottom 96-well plate. Take 5 $\mu$ l of each drug in column 1 of the original plate and transfer to column 1 in the master (round-bottom) plate. Next, take 5 $\mu$ l of each drug in column 2 of the original plate and transfer to column 1 of the master plate (same column as before). Repeat until 5 $\mu$ l from columns 1-6 of the original plate have been combined into column 1 in the master plate. Each master plate well will contain 30 $\mu$ l total.



Original plate

Master plate (30 $\mu$ l per well)

- Transfer the entire 30 $\mu$ l drug mix from each well on master plate to prepared 96-well plate of cells (containing 7000 cells in 100 $\mu$ l per well (for 5 $\mu$ M final concentration of each drug). Controls we used:

negative	positive	stock	dilute to:	
no treatment	5 $\mu$ M Celesterol	10mM	100 $\mu$ M	1 $\mu$ l in 99 $\mu$ l DME/10%FBS
DMSO (30 $\mu$ l)	35 $\mu$ M CdCl <sub>2</sub>	10mM	1mM	10 $\mu$ l in 90 $\mu$ l DME/10%FBS

- Incubate cells with drugs for 18hr.

#### Luciferase Assay

- Prepare Bright-Glo according to instructions. Aliquot into 15ml conical tubes (make ~10.5ml aliquots) and wrap in foil. Store aliquots at  $-70^{\circ}\text{C}$ .
- Remove cells from incubator and allow to sit at room-temp for 5min (to allow to cool to ~room temp as Bright-Glo reagent is temperature sensitive).
- Pour Bright-Glo reagent into multi-channel reservoir tray and add 100 $\mu$ l Bright-Glo reagent directly to each well containing 100  $\mu$ l of cells (Bright-Glo reagent is used at a volume of 1:1). *Note: be aware that the Bright-Glo luciferase reagent is time-sensitive so try to work as fast as possible and avoid having the reagent at room temp longer than necessary.*
- Incubate 5min at room temp to allow the cells to lyse (cover plate with foil as luciferase assay is light sensitive).
- Transfer 200 $\mu$ l of each well to white luminometer plate.
- Read plate immediately in the Lmax luminometer.
  - Turn on Lmax luminometer (switch is in the back of the instrument).
  - Insert the 96-well plate into the plate-slot.
  - Open the "Bright-Glo" plate reading program on the desktop.
  - Once the program is opened, click "READ". *Note: it takes about 15min for the luminometer to read a 96-well plate. Therefore, if you have multiple plates you should stagger the luciferase procedure so the reaction to be read by the plate reader immediately after the 5min incubation.*

#### Comments:

It is useful to keep in mind that pooling the compounds may hide the detection of a potentially positive response. If one drug out of the six in the combined pool is toxic enough to kill the cells (or prohibit a heat shock response) a very weak induction—if at all—will be detected.

We screen the drugs at 5 $\mu$ M, however, some drugs have been found to have wide-ranging optimization conditions.

**References:**

Matsumoto, Shuji et.al. Toxicology

Johnson, P.H. et. al, *Multiplex gene expression analysis for high-throughput drug discovery: screening and analysis of compounds affecting genes overexpressed in cancer cells*. Mol Cancer Ther. 2002 Dec;1(14):1293-304.