

RNAi FOR MAMMALIAN CELLS
(Silencer siRNA Construction Kit)

Materials/Reagents:

Silencer siRNA construction kit (Ambion cat #1620)
DNA oligos, sense and antisense. Design 4-6 sets for each target gene.
Oligofectamine transfection reagent (Invitrogen cat # 12252-011)
Cells, standard tissue culture supplies

Procedure:

Step 1: Design and order DNA oligos:

- Look for the first “AA” 75 nucleotides downstream of the transcription start site
- Include the following 19 additional nucleotides
- Verify that the GC content is between 30-70 percent (lower may be better)
- Perform a BLAST search to verify sequence specificity to target gene
- Determine sense and antisense DNA sequences and add 8 nt leader sequence “CCTGTCTC” to 3’ end according to the example below

(There is also a tool on the Ambion website that will do this for you)

Example:

Target cDNA sequence: 5’ AACGATTGACAGCAGCGGATTGCC 3’

Antisense siRNA template: 5’ AACGATTGACAGCGGATTGCCCTGTCTC 3’

Sense siRNA template: 5’ AAGGCAATCCGCTGTCAATCGCCTGTCTC 3’

Step 2: Resuspend oligos:

Resuspend oligos to approximately 200 μ M in nuclease-free dH₂O.

Dilute to 100 μ M stocks.

Step 3: Hybridization and fill-in

Mix (in separate tubes for sense and antisense):

Sense reaction

T7 promoter primer	2 μ l
DNA hybridization buffer	6 μ l
Sense template	<u>2 μl</u>
	10 μ l

Antisense reaction

T7 promoter primer	2 μ l
DNA hybridization buffer	6 μ l

To the ds RNA, add:

Digestion buffer	6 μ l
Nuclease-free dH ₂ O	48.5 μ l
RNase	3 μ l
DNase	2.5 μ l

Incubate 2 hours at 37°C.

Add 400 μ l siRNA binding buffer.

Incubate 2-5 min at room temp.

(Preheat 100 μ l/reaction dH₂O to 75°C)

Prewet filter cartridge with 100 μ l siRNA wash buffer.

Bind siRNA by adding to column and spinning at 10,000 rpm for 1 min.

Wash with 2X 500 μ l wash buffer.

Elute in 100 μ l 75°C dH₂O.

Quantitate siRNA and store at -80°C until ready to use.

Dilute to 20 μ M for siRNA transfection.

Step 6: siRNA transfection with Oligofectamine in 6 well plates

On the day before transfection, cells should be split into antibiotic-free media.

Cells should be 50% confluent on the day of transfection.

For each well of a 6 well plate:

Tube 1: Mix 12 μ l of 20 μ M siRNA with 200 μ l Opti-MEM (or serum free media)

Tube 2: Mix 12 μ l of Oligofectamine with 48 μ l Opti-MEM in a separate tube

Incubate 7 min at room temp

Combine tubes 1 and 2, mix by inversion (important), and incubate for 20 min at room temp.

Add 128 μ l fresh Opti-MEM for a final volume of 400 μ l.

Add 400 μ l of siRNA/Oligofectamine mix to the cells, swirl to mix.

After 48 hours, harvest cells to test for effectiveness of RNAi by RT-PCR or Western blot.

Comments:

For each gene, 4-6 target sites should be tested.

A nonspecific sequence can be used as a negative control.

Note for HSF1

I have used the target sequence aaatcaagaggaaagtgacc for HSF1 in MCF7 cells. Western analysis showed a 69 fold reduction in HSF1 protein. I still have to test for lack of HSF1 in a functional assay. So far in my hands, HeLa cells do not survive HSF1 RNAi.

References:

Silencer siRNA Construction Kit protocol (pdf file available online from Ambion)
siRNA user's guide: www.mpibpc.gwdg.de/abteilungen/100/105/

Submitted by: Sandy Westerheide