

LIGATION

This method works well for blunt end ligations, as well as sticky end ligations. Procedure from BRL:

1. Generally, 50 ng-100 ng of vector is used per ligation reaction.
2. Calculate from this amount of vector the amount of insert needed to achieve either a 1:1 or 2:1 molar ratio of insert:vector.
3. Add these amounts of vector and insert to Eppendorf tube in a combined volume of ≤ 15 μ l. (If volume exceeds this amount, vector and insert DNAs can be concentrated by ethanol precipitation.)
4. Add 4 μ l of 5x ligation buffer, made as recommended by BRL:
 - 250 mM Tris-HCl pH 7.6
 - 50 mM MgCl₂
 - 25% (w/v) polyethylene glycol 8000
 - 5 mM ATP
 - 5 mM DTT
5. Add 1-2 μ l T4 DNA ligase (BRL) for a final volume of 20 μ l.
 - * NOTE: 1 μ l is 1 unit of the BRL stock enzyme and this is plenty for 100 ng of vector and insert. Do not use 5 units of ligase as this is too much enzyme and the concentration of glycerol in the reaction will inhibit ligase activity.
6. Ligate at least 4 hours (to overnight) at room temperature.