

### LUCIFERASE ASSAY

1. Remove media from cells, wash once with 1xPBS, add 1 ml ice cold PBS to cells, scrape and transfer to eppendorf tube. Spin 15 sec at full speed.
2. Resuspend pellet in 150  $\mu$ l of 25 mM HEPES, pH 7.18, 15 mM MgSO<sub>4</sub>, 1 mM DTT.
3. Freeze and thaw 3 times using dry ice and a 37°C waterbath. After the final thaw, spin at full speed for 5 min, save supernatant. Do BioRad protein assay.
4. Take 30 to 50  $\mu$ g of each sample and add sample buffer to a volume of 175 $\mu$ l. Add 175  $\mu$ l of assay buffer. The total volume should be 350  $\mu$ l.

#### Sample Buffer

25 mM HEPES pH 7.8  
15 mM MgSO<sub>4</sub>

#### Assay Buffer

25 mM HEPES pH 7.8  
15 mM MgSO<sub>4</sub>  
5 mM ATP  
1  $\mu$ g/ml BSA  
(For 5 mls, add 50  $\mu$ l stock ATP and 5  $\mu$ l stock BSA)

5. Inject 100  $\mu$ l of luciferin into each tube using luminometer injection port and take a ten second integrated light reading. Be sure luciferin is at room temperature before using. Use 1  $\mu$ l of stock luciferase (1 ng) as a control. If everything is okay, a reading of 16,000,000 to 22,000,000 RLU should result.