

WESTERN BLOTTING ANALYSES USING ALKALINE PHOSPHATASE  
DETECTION

All steps should be done at room temperature with agitation.

- 1) Block filter for 1 hour in PBS + 5 mg/ml BSA -or- PBS + 2.5 % dry milk/0.2 % TWEEN-20.
- 2) Incubate with 1° antibody diluted in PBS/BSA -or- 1mg/ml BSA for 1 hour.
- 3) Wash 3 x 10 min with PBS + .2% TWEEN 20 (or 0.05% NP-40).
- 4) Incubate 30 min. - ~~Ab.~~ (with ~~2<sup>o</sup>~~ ~~anti~~ ~~goat~~ ~~antibody~~ ~~or~~ ~~anti-rabbit~~ Abs, dilutions should be 1:1000-1:2000 in PBS + 5 mg/ml BSA or 2.5 % dry milk.
- 5) Wash 3 x 10 min with PBS + 0.05% TWEEN (or NP-40).
- 6) \*\*Detect proteins by adding alkaline phosphatase buffer (100 mM tris pH 9.5; 100 mM NaCl; 5 mM MgCl<sub>2</sub>) containing substrate (for BioRad premixed substrates, use 6.6 µl NBT/ml and 3.3 µl BCIP/ml). HSP70 bands should show up within a couple minutes. Color reactions can be stopped by rinsing with dH<sub>2</sub>O.

\*\*Proteins can also be detected using peroxidase conjugated secondary antibodies. The drawback to this method is that the bands tend to fade over time. If using a peroxidase-conjugated 2°Ab instead of AP conjugated 2°: After post-incubation washes with 2°Ab, develop in 0.05% 45 chloro-1-naphthol, 0.01% H<sub>2</sub>O<sub>2</sub> in PBS. Once fully developed, stop with an H<sub>2</sub>O wash. Prepare solutions just before use: 10 ml PBS, 2 ml .3% chlororaphthol in MeOH, 4 µl 30% H<sub>2</sub>O<sub>2</sub>.