

PREPARATION OF DIDEOXY WORKING SOLUTIONS

(to sequence with radioactive dATP)

1. Prepare 0.5 mM Stocks of dCTP, dGTP, dTTP and 10 mM stocks of ddATP, ddTTP, ddCTP, ddGTP  
(These keep at -70°C indefinitely).
2. Prepare Zero Solutions.  
From 100 mM stock solutions, make 0.5 mM stocks (5 µl in 99.5 µl dH<sub>2</sub>O).

	A	C	G	T
0.5 mM dCTP	20 µl	1 µl	20 µl	20 µl
0.5 mM dGTP	20 µl	20 µl	1 µl	20 µl
0.5 mM dTT	20 µl	20 µl	20 µl	1 µl
1XTE	20 µl	20 µl	20 µl	20 µl

(These keep at -20°C about 2 months)

3. Prepare dideoxy nucleotide solutions

ddATP = 0.1 mM                      1 µl of 5 mM ddATP in 99 µl dH<sub>2</sub>O

ddCTP = 0.1 mM                      1 µl of 5 mM ddCTP in 99 µl dH<sub>2</sub>O

ddGTP = 0.3 mM                      6 µl of 5 mM ddGTP in 94 µl dH<sub>2</sub>O

ddTTP = 0.5 mM                      10 µl of 5 mM ddTTP in 90 µl dH<sub>2</sub>O

(These concentrations can be varied to give you optimal results for the sequence of interest)

(The solutions keep at -20°C about 2 months)

4. Combine ddNTP and dNTP  
Mix ddNTPs and dNTPs in a ratio of:

1:1 for longer products

2:1 for shorter products

(For sequencing, you will need 3 µl of each mix per clone sequenced.)

dsDNA System Nucleotide Mix formulations

component	C Nucleotide Mix	A Nucleotide Mix	T Nucleotide Mix	G Nucleotide Mix	Chase
ddCTP	66 $\mu$ M	---	---	---	---
ddATP	---	300 $\mu$ M	---	---	---
ddTTP	---	---	---117 $\mu$ M	---	--
ddGTP	---	---	---	66 $\mu$ M	---
dCTP	1.66 $\mu$ M	33 $\mu$ M	33 $\mu$ M	33 $\mu$ M	2 mM
dATP	---	---	---	---	2 mM
dTTP	33 $\mu$ M	33 $\mu$ M	1.66 $\mu$ M	33 $\mu$ M	2 mM
dGTP	33 $\mu$ M	33 $\mu$ M	33 $\mu$ M	1.66 $\mu$ M	2 mM
NaCl	50 mM	50 mM	50 mM	50 mM	50 mM
Tris-HCl pH 7.5	10 mM	10 mM	10 mM	10 mM	---
Tris-HCl pH 8.3	---	---	---	---	34 mM
MgCl <sub>2</sub>	10 mM	10 mM	10 mM	10 mM	6 mM
dithiothreitol	1 mM	1 mM	1 mM	1 mM	5 mM

*SEQUENCING SOLUTIONS*

Denaturing solution

500  $\mu$ l: 100  $\mu$ l 10 N NaOH

2  $\mu$ l 0.5 M EDTA  
398  $\mu$ l dH<sub>2</sub>O

Chase solution

1 ml: 20  $\mu$ l 100 mM dATP  
20  $\mu$ l 100 mM dCTP  
20  $\mu$ l 100 mM dGTP  
20  $\mu$ l 100 mM dTTP  
920  $\mu$ l dH<sub>2</sub>O

10x Hybridization buffer

1 ml: 200  $\mu$ l 1 M Tris pH 7.5  
2  $\mu$ l 0.5 M EDTA  
125  $\mu$ l 4 M NaCl  
100  $\mu$ l 1 M MgCl<sub>2</sub>  
20  $\mu$ l 0.5 M DTT  
553  $\mu$ l dH<sub>2</sub>O

6x proteinase K buffer

500 ml: 30  $\mu$ l 1 M Tris pH 7.6  
60  $\mu$ l 0.5 M EDTA  
25.43 g LiCl  
30 ml 20% SDS

SEALING GEL

10 ml: 8 ml dH<sub>2</sub>O  
2 ml 40% Acrylamide  
--10 mg APS  
--30  $\mu$ l TEMED

SEQUENCING GEL

	<u>8% Acrylamide gel</u>	<u>6% Acrylamide gel</u>
50 ml:	10 ml 40% Acrylamide	7.5 ml 40% Acrylamide
	10 ml dH <sub>2</sub> O	12.5 ml dH <sub>2</sub> O
	25 g urea	25 g urea
	10 ml 5xTBE	10 ml 5xTBE
	--50 mg APS	--50 mg APS
	--22 $\mu$ l TEMED	--22 $\mu$ l TEMED

