

SITE DIRECTED MUTAGENESIS

This technique can be used to generate single point mutations or to delete/insert single or multiple amino acids.

Procedure:

1. PCR reaction
5 μ l 10X PFU buffer
X μ l (5-50ng) dsDNA template
X μ l (125ng) primer #1
X μ l (125ng) primer #2
X μ l 10X dNTPs
ddH₂O to a final volume of 50 μ l
1 μ l PFU polymerase

Cycling parameters:

Segment	Cycles	Temperature	Time
1	1	95°C	30 seconds
2	12-18	95°C	30 seconds
		55°C	1 minute
		68°C	2minutes/kb of plasmid length

Type of mutation desired	Number of cycles
Point mutations	12
Single amino acid changes	16
Multiple amino acid changes	18

2. Chill PCR reaction on ice and add 1 μ l DpnI. Mix and incubate at 37°C for 1 hour to digest the parental supercoiled dsDNA.
3. Transform 1 μ l of the DpnI treated DNA and plate on appropriate medium. Screen colonies for desired mutants.

Comments:

1. It is very important that the primers be designed to have a T_m=78.
2. It is useful to set up a series of PCR reactions with 5,10,20, and 50ng of template DNA to determine optimal concentration.

Reference:

QuikChange Site Directed Mutagenesis Kit Manual at www.stratagene.com