

DUAL LUCIFERASE REPORTER SYSTEM (PROMEGA)  
(24 well plate)

Transfections

Cells should be transfected with a test reporter (firefly luciferase) and a normalization reporter driven by a constitutive promoter (renilla luciferase)  
PRL-TK (thymidine kinase promoter) works well at a ratio of 1:200.

Cell Lysis

1. Thaw 5X PLB and prepare appropriate amount of 1X PLB by diluting in dH<sub>2</sub>O.
2. Remove media from cells by aspiration.
3. GENTLY wash with 1X PBS.
4. Add 100  $\mu$ l 1X PLB to each well. Place on shaker at RT for 30 min.

Reagent Preparation

1. Prepare Luciferase assay reagent II:

Resuspend lyophilized Luciferase Assay Substrate in 10 mL of Luciferase Assay Buffer II and label "LARII" (stable for one year at -80).

\*need 100  $\mu$ l per assay, plus about 2 mL to tube volume for luminometer

2. Prepare Stop and Glo:

a) Make 50X Stop and Glo substrate: Transfer 200  $\mu$ l Stop and Glo substrate solvent into lyophilized Stop and Glo substrate, vortex (stable for one year at -80).

b) Make 1X Stop and Glo substrate: Add 1 vol of 50X Stop and Glo substrate to 50 vol Stop and Glo buffer. Store in dark glass bottle.

\* need 100  $\mu$ l per assay, plus about 2 mLs to tube volume for luminometer

Luminescence measurement

Read samples in Molecular Devices "LMax" luminometer:

1. Dispense 20  $\mu$ L cell lysate into white opaque luminometer microtiter plate.
2. Prime P and M injectors (7 injections each) with LARII and Stop and Glo.
3. Inject 100  $\mu$ l LARII per well with P-injector. Read with a 2 sec premeasurement delay, 10 sec integration to detect firefly luciferase.
4. Inject 100  $\mu$ l Stop and Glo per well with the M injector. Read with a 2 sec premeasurement delay, 10 sec integration to detect renilla luciferase. Print out data.

Wash injectors

1. Place both tubes into a 50ml Falcon tube containing ddH<sub>2</sub>O. Select **Wash Injectors...** from the **Control** menu. Accept the default settings (30 injections) and click **OK**.
2. Repeat this wash step with 70%EtOH (30 injections). Replace each injector tube into the empty Falcon tubes and perform a final “Wash” that will just dry the injector tubing.

**Submitted by:** Sandy Westerheide