

PREPARATION OF COMPETENT *E. coli* CELLS FOR TRANSFORMATION

Shawn's Protocol

*Large scale stocks of *E. coli* cells can be made competent for transformation with plasmid DNAs, and then stored frozen at -70°C.*

- 1) Prepare a streak-plate of DH-1 cells on an LB plate. Incubate at 37°C overnight.

Note: This protocol is for 1-liter of competent DH-1 cells. Will yield 50 ml competent cells to be aliquotted into 200 µl or smaller aliquots.

Have prepared: 4 x 250 ml L-broth in 1-liter flasks
 1100 ml ice-cold 0.1 M CaCl₂

- 2) Inoculate two 5 ml L-broth cultures (in 15 ml snap-cap tubes) with a single DH-1 colony and shake overnight at 37°C in a New Brunswick G-25 shaker at 400 rpm.
- 3) Transfer the overnight cultures equally into four 1-liter flasks containing 250 ml L-broth each. Shake the flasks at 37°C until the bacteria reaches OD₅₅₀ = 0.5, which will normally require approximately 2 hours.
- 4) Cool the bacteria on ice to roughly 5°C, then transfer to four 500 ml centrifuge bottles, and centrifuge in a Sorvall GSA rotor at 4°C for 10 minutes at 5000 rpm (10,400 x g).
- 5) Resuspend bacterial pellets in 200 ml ice-cold 0.1 M CaCl₂. Combine cells into a single bottle, and centrifuge again as above.
- 6) Resuspend pellet in 250 ml ice-cold 0.1 M CaCl₂, and keep on ice for 4-8 hours.
Note: Competent cells can be obtained with incubation times on ice as low as 20 minutes, but the longer times greatly increase the overall transformation efficiency.
- 7) Recentrifuge cells as above. Resuspend pellet in 43 ml ice-cold 0.1 M CaCl₂ mixed with 7 ml sterile glycerol.
- 8) Dispense cells into 200 µl aliquots in 1.5 ml eppendorf tubes. *Aliquot in the cold room to prevent the cells from warming.*
- 9) Store the cells frozen at -70°C until use. Bacteria prepared in this manner can be successfully used after up to one year's storage and will still maintain a high transformation efficiency.