

## BIORAD ASSAY PROTEIN CONCENTRATION DETERMINATION

Bio-Rad Protein Assay Dye Reagent Concentrate is a colorimetric assay for protein concentration. Similar to the Lowry assay, but with the following improvements: the reaction reaches 90% of its maximum color development within 15 minutes and the color changes not more than 5% in 1 hour. The assay is based on the reaction of proteins with an alkaline copper tartrate solution and Folin reagent. A standard curve should be prepared each time the assay is performed.

### **Materials:**

Bio-Rad Protein Assay Dye Reagent Concentrate (cat#500-0006)

BSA standard (2 $\mu$ g/ $\mu$ l)

Disposable cuvettes (BioRad #223-9955)

Spectrophotometer (we use a Hitachi U-2000)

### **Procedure:**

- 1) Dilute the BioRad assay dye 1:5 and mix well. Make enough for 1ml/sample plus 7ml for the standards. (note: it is useful to make enough for 2ml for each unknown sample plus a few extra mls in case any samples need to re-checked).

10ml

2 ml BioRad dye

8 ml dH<sub>2</sub>O

- 2) Prepare standards (for best results, the standards should always be prepared in the same buffer as the sample).
  - Dispense 1 ml aliquots of the diluted dye into 7 tubes.
  - Add the appropriate amount of BSA standard to each tube: (0, 0, 2, 4, 6, 8, 10 $\mu$ g)
  - Mix each tube well; there should be a gradient of blue color ranging from the lowest to highest amount of BSA added.
- 3) Prepare sample:
  - Dispense 1ml of the diluted dye to the number of tubes needed for each sample.
  - For each sample, add 1 $\mu$ l and 2 $\mu$ l respectively to 2 separate tubes.
  - Mix each sample well.
- 4) Transfer each sample to a disposable plastic cuvette.
- 5) Read each sample in the spectrophotometer (595l):
  1. Press the "Main Menu" button on the spec.
  2. Select #6 (Test Menu). ENTER
  3. Select #1 (Load Test). ENTER.
  4. Select #15 (BioRad AS). ENTER.
  5. Read the BSA Standards first.

- Place one of the “0” cuvettes in the “blank” chamber and place the other “0” cuvette in the “sample” chamber. Press “Auto Zero”. To begin reading, press “Start”.
  - Remove the “0” cuvette from the sample chamber and place the “2  $\mu\text{g}$ ” sample there. Press “Start”.
  - Repeat this with the 4, 6, 8, and 10 $\mu\text{g}$  sample.
  - A graph will immediately print out following the 10 $\mu\text{g}$  sample. If your standard curve should be 0.995 or better to be considered accurate.
6. Place the first unknown sample in the sample chamber. Press “start”. The concentration will immediately appear; the unit of the sample will be \_\_\_ $\mu\text{g}$ / how many  $\mu\text{l}$  of unknown sample you used.